

New  
Innovation

# Drummond

## Portable Pipet-Aid® XL

### New Ergonomic Design Reduces Fatigue When Pipetting Under a Hood

- Longer lightweight handle lowers your arm position to reduce strain
- New quiet, more efficient pump extends operation period between charges
- Can be charged while in use to eliminate downtime
- Three-speed settings provide more precise control for both aspirate and dispense modes

Pipetting under a hood doesn't need to be a pain in the neck. The new Drummond Portable Pipet-Aid XL has a longer handle to reduce the amount of arm lift required to perform the same operation as with a conventional pipettor. The ergonomic design includes an adjustable sliding hand rest and a removable stand which enables the unit to be set down without contaminating the pipet. The Portable Pipet-Aid XL also features a new quieter, more powerful pump, and a new EZ Grip™ nosepiece. The unit is supplied complete with power supply charger, 4 extra filters, and a holster-type wall bracket.



### DRUMMOND PORTABLE PIPET-AID® XL

| Cat. No.  | Description                |
|-----------|----------------------------|
| 4-000-105 | Portable Pipet-Aid XL 110V |

#### Replacement Parts

| Cat. No.  | Description                   |
|-----------|-------------------------------|
| 4-000-002 | Rubber Insert (2/pk)          |
| 4-000-051 | TC Filters (5/pk)             |
| 4-000-057 | TC Gaskets (pair)             |
| 4-000-036 | Battery Pack (XL only)        |
| 4-000-065 | XL Power Supply/Charger, 110V |
| 4-000-079 | EZ Grip Nosepiece             |

Compared to a conventional pipettor, the Portable Pipet-Aid XL lowers hand height required for the same task



*The Developers of the Original Pipet-Aid*

For more information, call 800/523-7480  
or visit our web site at [www.drummondsci.com](http://www.drummondsci.com)



**DRUMMOND**  
SCIENTIFIC COMPANY

500 Parkway, Box 700  
Broomall, PA 19008  
Call Toll Free: 1-800/523-7480  
1-610/353-0200  
FAX: 1-610/353-6204  
Web site: [www.drummondsci.com](http://www.drummondsci.com)