

NEW PRODUCT INNOVATIONS FROM DRUMMOND



**Look...
One
Hand!**



DECAPITATOR 2™*

Fast, Easy Microcentrifuge Tube Cap Remover

The new Drummond DECAPITATOR 2 provides a simple, secure one-handed method for uncapping microcentrifuge tubes. Simply slide the cap into the appropriate slot and apply a slight downward pressure. The DECAPITATOR 2 can accommodate domed cap PCR tubes and most tube sizes from 0.2 mL to 2.2 mL. Comes complete with a compact, sturdy lab bench stand.

DRUMMOND DECAPITATOR 2

Catalog Number	Quantity
6-000-002	1

*Patent 6,307,603 B1

Available through your local distributor



**DRUMMOND
SCIENTIFIC COMPANY**

500 Parkway, Box 700
Broomall, PA 19008
Call Toll Free: 1-800/523-7480
1-610/353-0200
FAX: 1-610/353-6204
Website: www.drummondsci.com



MADE IN USA

©Copyright, Drummond Scientific Company, 1999 Rev 11/01 Printed in USA
DECAPITATOR 2™ is a trademark of Drummond Scientific Company.