

# NEW PRODUCT INNOVATIONS FROM DRUMMOND



Mylar wrapped HEMATO-CLAD tubes help prevent contamination and provide an added layer of protection.

## HEMATO-CLAD™ Mylar® Wrapped Hematocrit Tubes

**Complies with FDA, NIOSH, CDC,  
and OSHA Safety Recommendations  
for the Use of Glass Capillary Tubes  
for Blood Collection\*** February 22, 1999

Drummond HEMATO-CLAD Mylar wrapped hematocrit tubes meet government agency safety recommendations for the use of glass capillary tubes for the collection of blood. HEMATO-CLAD tubes combine the precision of glass with the safety provided by the Mylar overwrap.

With HEMATO-CLAD tubes the chance of contamination due to glass fragments or aerosols is greatly reduced. Accidental breakage is significantly reduced particularly when pressing the tube into sealing clay. In the event of breakage, the Mylar wrapping will help contain the sample (blood) and provides an added layer of protection against cuts from broken glass.

**Precision and Safety...  
The Exclusive HEMATO-CLAD Combination**

### DRUMMOND HEMATO-CLAD TUBES

Catalog Number	Description	Package Size
1-000-7500-C	Plain (Blue)	200/Vial
1-000-7500-HC	Heparin Coated (Red)	200/Vial
1-000-750-PCC	Precalibrated Tube	100/Vial

**Available through your local distributor**



MADE IN USA

HEMATO-CLAD™ is a trademark of Drummond Scientific Company.  
Mylar® is a registered trademark of DuPont.  
©Copyright, Drummond Scientific Company, 1999  
Printed in USA Rev. 3/99

\* [www.fda.gov/cdrh/safety.html](http://www.fda.gov/cdrh/safety.html)  
[www.cdc.gov/niosh/homepage.html](http://www.cdc.gov/niosh/homepage.html)  
[www.osha-slc.gov/sltc/needlestick/](http://www.osha-slc.gov/sltc/needlestick/)



**DRUMMOND  
SCIENTIFIC COMPANY**

500 Parkway, Box 700  
Broomall, PA 19008  
Call Toll Free: 1-800/523-7480  
1-610/353-0200  
FAX: 1-610/353-6204  
Website: [www.drummondsci.com](http://www.drummondsci.com)